

The logo of Galgotias University is a stylized 'G' composed of several curved, overlapping bands in shades of yellow, orange, and blue. The word 'Mounting' is written in a large, bold, black sans-serif font, centered over the logo.

Mounting

GALGOTIAS
UNIVERSITY

Mounting




GALGOTIAS
UNIVERSITY

Mounting

- **Two types of mounting media**
 - Water based mounting media
 - Resinous mounting media
- **DPX & Canada balsam commonly used**

GALGOTIAS
UNIVERSITY

Mounting



- **To mount a slide**

- A. Apply drops of mounting medium upon tissue section.
- B. Hold coverslip at 45° allowing the drop to spread along the edge of the slip.
- C. Let go of slip and allow the medium to spread slowly.

- **Allow it to dry and section viewed under microscope**

Mounting



Mounting



Reference

- <https://www.leicabiosystems.com/knowledge-pathway/an-introduction-to-specimen-processing/>
- <https://tissuesampling.weebly.com/processing.html>

GALGOTIAS
UNIVERSITY